

## Product Datasheet

### Human IL-17RA Protein (Biotin) (orb1516405)

<b>Description</b>	Recombinant Human IL-17RA Protein conjugated with Biotin
<b>Endotoxins</b>	< 1.0 EU/μg as determined by LAL
<b>Target</b>	IL17RA
<b>Concentration</b>	> 0.5 mg/ml
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	C-His-Avi
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Protein Sequence</b>	33-320/866
<b>Purity</b>	> 95% as determined by Tris-Bis PAGE; > 92% as determined by SEC-HPLC
<b>MW</b>	Due to glycosylation, the protein migrates to 60-70 kDa based on Tris-Bis PAGE result. kDa
<b>Uniprot ID</b>	<a href="#">Q96F46</a>
<b>Expiration Date</b>	6 months from date of receipt.