

## Product Datasheet

### Human FGF21 Protein (orb1516454)

<b>Description</b>	Recombinant Human FGF21 Protein
<b>Endotoxins</b>	< 1.0 EU/μg as determined by LAL
<b>Target</b>	FGF21
<b>Form/Appearance</b>	Lyophilized powder
<b>Concentration</b>	> 0.5 mg/ml
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	N-mFc-Avi
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Protein Sequence</b>	29-209/209
<b>Purity</b>	> 95% as determined by Tris-Bis PAGE; > 95% as determined by SEC-HPLC
<b>MW</b>	Due to glycosylation, the protein migrates to 52-60 kDa based on Tris-Bis PAGE result. kDa
<b>Uniprot ID</b>	<a href="#">Q9NSA1</a>
<b>Expiration Date</b>	6 months from date of receipt.