

Product Datasheet

Human DLL4 Protein (orb1516466)

| | |
|--------------------------|--|
| Description | Recombinant Human DLL4 Protein |
| Endotoxins | < 1.0 EU/μg as determined by LAL |
| Target | DLL4 |
| Form/Appearance | Lyophilized powder |
| Concentration | > 0.5 mg/ml |
| Storage | Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles. |
| Tag | C-His-Avi |
| Note | For research use only |
| Application notes | Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water. |
| Protein Sequence | 27-524/685 |
| Purity | > 95% as determined by Tris-Bis PAGE; > 95% as determined by SEC-HPLC |
| MW | Due to glycosylation, the protein migrates to 62-68 kDa based on Tris-Bis PAGE result. kDa |
| Uniprot ID | Q9NR61 |
| Expiration Date | 6 months from date of receipt. |