

Product Datasheet

Human GM-CSF R Alpha Protein (Biotin) (orb1516469)

Description	Recombinant Human GM-CSF R Alpha Protein conjugated with Biotin
Endotoxins	< 1.0 EU/μg as determined by LAL
Target	CSF2RA
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His-Avi
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	23-320/400
Purity	> 95% as determined by Tris-Bis PAGE; > 95% as determined by SEC-HPLC
MW	Due to glycosylation, the protein migrates to 68-75 kDa based on Tris-Bis PAGE result. kDa
Uniprot ID	P15509
Expiration Date	6 months from date of receipt.