

## Product Datasheet

# Human /Monkey/rhesus macaque CD28 Protein (Biotin) (orb1516479)

<b>Description</b>	Recombinant Human /Monkey/rhesus macaque CD28 Protein conjugated with Biotin
<b>Endotoxins</b>	< 1.0 EU/μg as determined by LAL
<b>Target</b>	CD28
<b>Concentration</b>	> 0.5 mg/ml
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	N-His-Avi
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Protein Sequence</b>	17-152/220
<b>Purity</b>	> 95% as determined by Tris-Bis PAGE; > 95% as determined by SEC-HPLC
<b>MW</b>	Due to glycosylation, the protein migrates to 38-50 kDa based on Tris-Bis PAGE result. kDa
<b>Uniprot ID</b>	<b>P31041</b>
<b>Expiration Date</b>	6 months from date of receipt.