

Product Datasheet

Human Angiopoietin-2 Protein (Biotin) (orb1516535)

| | |
|--------------------------|---|
| Description | Recombinant Human Angiopoietin-2 Protein conjugated with Biotin |
| Endotoxins | < 1.0 EU/μg as determined by LAL |
| Target | ANGPT2 |
| Concentration | > 0.5 mg/ml |
| Storage | Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles. |
| Tag | N-His-Avi |
| Note | For research use only |
| Application notes | Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water. |
| Protein Sequence | 275-496/496 |
| Purity | > 95% as determined by Tris-Bis PAGE; > 95% as determined by SEC-HPLC |
| MW | 35 kDa |
| Uniprot ID | O15123 |
| Expiration Date | 6 months from date of receipt. |