

Product Datasheet

Human IL-8 Protein (orb1516542)

Description	Recombinant Human IL-8 Protein
Endotoxins	Not analyzed
Target	CXCL8
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Shipped at 4°C. Store at -20°C or -70°C for one year. Avoid repeated freeze/thaw cycles.
Tag	C-Fc
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	23-99/99
Purity	> 66% as determined by SDS-PAGE
MW	50-55 kDa
Uniprot ID	P10145
Expiration Date	6 months from date of receipt.