

Product Datasheet

Recombinant mouse IgG1 Fc protein (Isotype Control) (HEK293) (orb1516543)

Description	Recombinant mouse IgG1 Fc protein (Isotype Control) (HEK293)
Endotoxins	< 1.0 EU/μg as determined by LAL
Target	Ig gamma-1 chain C region secreted form
Preservatives	PBS (pH7.4).
Form/Appearance	Liquid
Concentration	□0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	Tag Free
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	98-324/324
Clonality	Recombinant
Purity	> 95% as determined by SDS-PAGE
MW	32 kDa
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)