

Product Datasheet

Mouse IgG1 Fc Protein (Isotype Control) (HEK293) (orb1516543)

Description	Recombinant Mouse IgG1 Fc Protein (Isotype Control) (HEK293)
Endotoxins	< 1.0 EU/μg as determined by LAL
Target	IGHG1
Form/Appearance	Liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	Tag Free
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	98-324/324
Clonality	Recombinant
Purity	> 95% as determined by SDS-PAGE
MW	32 kDa
Uniprot ID	P01868
Expiration Date	6 months from date of receipt.