

Product Datasheet

Human CD90 (THY1) Protein (orb1516548)

Description	Recombinant Human CD90 (THY1) Protein
Endotoxins	Not analyzed
Target	THY1
Form/Appearance	Liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	1-130/161
Purity	> 90% as determined by SDS-PAGE
MW	22-30 kDa
Uniprot ID	P04216
Expiration Date	6 months from date of receipt.