

Product Datasheet

Human CD4 Protein (orb1516549)

Description	Recombinant Human CD4 Protein
Endotoxins	Not analyzed
Target	CD4
Preservatives	Supplied as solution form in PBS or lyophilized from PBS (pH=7.4).
Form/Appearance	Liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	1-390/458
Purity	> 90% as determined by SDS-PAGE
MW	48 kDa
Uniprot ID	P01730
Expiration Date	6 months from date of receipt.