

Product Datasheet

Human ACE2 Protein (orb1516559)

Description	Recombinant Human ACE2 Protein
Endotoxins	Not analyzed
Target	ACE2
Preservatives	Lyophilized from 0.22 um filtered solution in PBS (pH=7.4). Normally 5% trehalose is added as protectant before Lyophilized.
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-hFc
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	18-740/805
Purity	> 95% as determined by SDS-PAGE
MW	115-130 kDa
Uniprot ID	Q9BYF1
Expiration Date	6 months from date of receipt.