

Product Datasheet

SARS-Cov-2 Spike trimer Protein (Biotin) (orb1516561)

Description	Recombinant SARS-Cov-2 Spike trimer Protein conjugated with Biotin
Endotoxins	< 1.0 EU/μg as determined by LAL
Target	SARS-CoV-2
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His-Avi
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	16-1188/1273
Purity	> 95% as determined by Tris-Bis PAGE
MW	Due to glycosylation, the protein migrates to 175-230 kDa based on Tris-Bis PAGE result. kDa
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.