

## Product Datasheet

# SARS-Cov-2 Spike RBD Protein (Biotin) (orb1516574)

<b>Description</b>	Recombinant SARS-Cov-2 Spike RBD Protein conjugated with Biotin
<b>Endotoxins</b>	Not analyzed
<b>Target</b>	SARS-CoV-2
<b>Concentration</b>	> 0.5 mg/ml
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	His & Avi
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Protein Sequence</b>	319-532/1273
<b>Purity</b>	> 95% as determined by SDS-PAGE
<b>MW</b>	34 kDa
<b>Uniprot ID</b>	<a href="#">P0DTC2</a>
<b>Expiration Date</b>	6 months from date of receipt.