

Product Datasheet

Human ACE2 Protein (Biotin) (orb1516577)

Description	Recombinant Human ACE2 Protein conjugated with Biotin
Endotoxins	Not analyzed
Target	SARS-CoV-2
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	His & Avi
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	18-740/805
Purity	> 95% as determined by SDS-PAGE
MW	100-110 kDa
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.