

Product Datasheet

MPXV B6R Protein (orb1516580)

Description	Recombinant MPXV B6R Protein
Endotoxins	Not analyzed
Target	Monkeypox virus
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	64-276/317
Purity	> 90% as determined by SDS-PAGE
MW	25.7 kDa
Uniprot ID	Q8V502
Expiration Date	6 months from date of receipt.