

## Product Datasheet

### MPXV B14R Protein (orb1516582)

<b>Description</b>	Recombinant MPXV B14R Protein
<b>Endotoxins</b>	Not analyzed
<b>Target</b>	Monkeypox virus
<b>Form/Appearance</b>	Lyophilized or liquid
<b>Concentration</b>	> 0.5 mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Protein Sequence</b>	21-303/326
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>MW</b>	34.3 kDa
<b>Uniprot ID</b>	<a href="#">Q8V502</a>
<b>Expiration Date</b>	6 months from date of receipt.