

## Product Datasheet

### Human MCP-1 Protein (orb1516587)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Recombinant Human MCP-1 Protein  |
| <b>Endotoxins</b>        | Not analyzed   |
| <b>Target</b>            | CCL2   |
| <b>Preservatives</b>     | PBS (pH=7.4)   |
| <b>Form/Appearance</b>   | Lyophilized or liquid  |
| <b>Concentration</b>     | > 0.5 mg/ml  |
| <b>Storage</b>           | Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.   |
| <b>Tag</b>               | C-His-Avi  |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | Centrifuge tubes before opening.<br>Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water. |
| <b>Protein Sequence</b>  | 24-99/148  |
| <b>Purity</b>            | > 90% as determined by SDS-PAGE  |
| <b>MW</b>                | 15/18 kDa  |
| <b>Uniprot ID</b>        | <a href="#">P13500</a>   |
| <b>Expiration Date</b>   | 6 months from date of receipt.   |