

Product Datasheet

Human TGFBR2 Protein (orb1516589)

Description	Recombinant Human TGFBR2 Protein
Endotoxins	Not analyzed
Target	Tgfr2
Preservatives	PBS (pH=7.4)
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	hFc & His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	23-159/567
Purity	> 90% as determined by SDS-PAGE
MW	50-70 kDa
Uniprot ID	P37173
Expiration Date	6 months from date of receipt.