

Product Datasheet

Mouse IL-4 Protein (orb1516594)

Description	Recombinant Mouse IL-4 Protein
Endotoxins	Not analyzed
Target	IL4
Preservatives	PBS (pH=7.4)
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	1-140/140
Purity	> 95% as determined by SDS-PAGE
MW	21 kDa
Uniprot ID	P07750
Expiration Date	6 months from date of receipt.