

Product Datasheet

Mouse IL-21 Protein (orb1516602)

Description	Recombinant Mouse IL-21 Protein
Endotoxins	Not analyzed
Target	IL21
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His Tag
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	18-146/146
Purity	> 90% as determined by SDS-PAGE
MW	16.8 kDa
Uniprot ID	Q9HBE4
Expiration Date	6 months from date of receipt.