

Product Datasheet

Human Adiponectin Protein (orb1516609)

Description	Recombinant Human Adiponectin Protein
Endotoxins	Not analyzed
Target	ADIPOQ
Form/Appearance	Lyophilized or liquid
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His Tag
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	19-244/244
Purity	> 90% as determined by SDS-PAGE
MW	26.5 kDa
Uniprot ID	Q60994
Expiration Date	6 months from date of receipt.