

Product Datasheet

Human NT-proBNP Protein (orb1516611)

Description	Recombinant Human NT-proBNP Protein
Endotoxins	Not analyzed
Target	NPPB
Form/Appearance	Lyophilized or liquid
Concentration	> 1.0 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	N-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	27-134/134
Purity	> 90% as determined by SDS-PAGE
MW	14 kDa
Uniprot ID	P16860
Expiration Date	6 months from date of receipt.