

Product Datasheet

Human IL-6R-Beta Protein (orb1516616)

Description	Recombinant Human IL-6R-Beta Protein
Endotoxins	Not analyzed
Target	IL6ST
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-hFc-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	23-618/618
Purity	> 90% as determined by SDS-PAGE
MW	95.2 kDa
Uniprot ID	P40189
Expiration Date	6 months from date of receipt.