

Product Datasheet

Human LH Alpha Protein (orb1516620)

Description	Recombinant Human LH Alpha Protein
Endotoxins	Not analyzed
Target	LHB
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	25-116/116 (Alpha chain) & 21-141/141 (Beta chain)
Purity	> 90% as determined by SDS-PAGE
Uniprot ID	P01229
Expiration Date	6 months from date of receipt.