

Product Datasheet

Human TSH Alpha/Beta Heterodimer Protein (orb1516621)

Description	Recombinant Human TSH Alpha/Beta Heterodimer Protein
Endotoxins	Not analyzed
Target	TSHB
Form/Appearance	Lyophilized powder
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	25-116/116 (Alpha chain) & 21-138/138 (Beta chain)
Purity	> 90% as determined by SDS-PAGE
Uniprot ID	P01222 , P01215
Expiration Date	6 months from date of receipt.