

Product Datasheet

SARS-Cov-2 Protein (orb1516626)

Description	Recombinant SARS-Cov-2 Protein
Endotoxins	Not analyzed
Target	SARS-CoV-2
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	16-1213/1273
Purity	> 90% as determined by SDS-PAGE
MW	150-180 kDa
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.