

Product Datasheet

SARS-Cov-1 Spike S1 Protein (orb1516634)

Description	Recombinant SARS-Cov-1 Spike S1 Protein
Endotoxins	Not analyzed
Target	SARS
Preservatives	PBS (pH=7.4)
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	14-667/1255
Purity	> 90% as determined by SDS-PAGE
MW	110-140 kDa
Uniprot ID	P59594
Expiration Date	6 months from date of receipt.