

Product Datasheet

Human IL-7R Alpha Protein (orb1516641)

Description	Recombinant Human IL-7R Alpha Protein
Endotoxins	Not analyzed
Target	IL7R
Preservatives	PBS (pH=7.4)
Form/Appearance	Liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	21-239/459
Purity	> 95% as determined by SDS-PAGE
MW	40-70 kDa
Uniprot ID	P16871
Expiration Date	6 months from date of receipt.