

Product Datasheet

SARS-Cov-2 Spike S1 Protein (orb1516642)

Description	Recombinant SARS-Cov-2 Spike S1 Protein
Endotoxins	Not analyzed
Target	SARS-CoV-2
Preservatives	PBS (pH=7.4)
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles
Tag	His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	14-685/1273
Purity	> 90% as determined by SDS-PAGE
MW	120-140 kDa
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.