

Product Datasheet

Mouse CD40L Protein (orb1516651)

Description	Recombinant Mouse CD40L Protein
Endotoxins	Not analyzed
Target	CD40LG
Preservatives	PBS (pH=7.4)
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Tag	His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	112-260/260
Purity	> 95% as determined by SDS-PAGE
MW	20 kDa
Uniprot ID	P27548
Expiration Date	6 months from date of receipt.