

Product Datasheet

Human TNFSF13B Protein (orb1516652)

Description	Recombinant Human TNFSF13B Protein
Endotoxins	Not analyzed
Target	TNFSF13B
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Tag	N-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	134-285/285
Purity	> 95% as determined by SDS-PAGE
MW	21 kDa
Uniprot ID	Q9Y275
Expiration Date	6 months from date of receipt.