

Product Datasheet

Human NMP-22 Protein (orb1516653)

Description	Recombinant Human NMP-22 Protein
Endotoxins	Not analyzed
Target	NUMA1
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	N-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	2-410/2115
Purity	> 95% as determined by SDS-PAGE
MW	52 kDa
Uniprot ID	Q14980
Expiration Date	6 months from date of receipt.