

Product Datasheet

Human Apolipoprotein E Protein (orb1516656)

Description	Recombinant Human Apolipoprotein E Protein
Endotoxins	Not analyzed
Target	APOE
Preservatives	PBS (pH=7.4)
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	19-317/317
Purity	> 95% as determined by SDS-PAGE
MW	40 kDa
Uniprot ID	P02649
Expiration Date	6 months from date of receipt.