

Product Datasheet

Human MMP-3 Protein (orb1516657)

Description	Recombinant Human MMP-3 Protein
Endotoxins	Not analyzed
Target	MMP3
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	18-477/477
Purity	> 95% as determined by SDS-PAGE
MW	60 kDa
Uniprot ID	P08254
Expiration Date	6 months from date of receipt.