

Product Datasheet

Human IL-6R Alpha Protein (Bio) (orb1516658)

Description	Recombinant Human IL-6R Alpha Protein conjugated with Bio
Endotoxins	Not analyzed
Target	IL6R
Preservatives	PBS (pH=7.4)
Form/Appearance	Liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	1-468/468
Purity	> 95% as determined by SDS-PAGE
Uniprot ID	P08887
Expiration Date	6 months from date of receipt.