

Product Datasheet

Human IGFBP1 Protein (orb1516666)

| | |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Recombinant Human IGFBP1 Protein |
| Endotoxins | Not analyzed |
| Target | IGFBP1 |
| Form/Appearance | Lyophilized or liquid |
| Concentration | > 0.5 mg/ml |
| Storage | Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles. |
| Tag | C-His |
| Note | For research use only |
| Application notes | Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water. |
| Protein Sequence | 26-259/259 |
| Purity | > 95% as determined by SDS-PAGE |
| MW | 34 kDa |
| Uniprot ID | P08833 |
| Expiration Date | 6 months from date of receipt. |