

Product Datasheet

Mouse Thy-1 Protein (orb1516669)

Description	Recombinant Mouse Thy-1 Protein
Endotoxins	Not analyzed
Target	THY1
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	20-131/162
Purity	> 95% as determined by SDS-PAGE
MW	23-30 kDa
Uniprot ID	P04216
Expiration Date	6 months from date of receipt.