

## Product Datasheet

### SARS-Cov-2 N Protein (orb1516693)

<b>Description</b>	Recombinant SARS-Cov-2 N Protein
<b>Endotoxins</b>	Not analyzed
<b>Target</b>	SARS-CoV-2
<b>Preservatives</b>	20 mM Tris-HCl (pH=8.0)
<b>Form/Appearance</b>	Lyophilized or liquid
<b>Concentration</b>	> 0.5 mg/ml
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	His
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Protein Sequence</b>	1-419/419
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>MW</b>	50 kDa
<b>Uniprot ID</b>	<a href="#">P0DTC2</a>
<b>Expiration Date</b>	6 months from date of receipt.