

Product Datasheet

Recombinant human IL-20R alpha protein, N-His (orb1516743)

Catalog Number	orb1516743
Category	Proteins
Description	Recombinant human IL-20R alpha protein, N-His, with a concentration of >0.5 mg/ml, expressed in E.coli, with a molecular weight of 28 kDa, and a purity of >95% by SDS-PAGE. It was purified by AC, and has a N-His tag.
Target	Interleukin-20 receptor subunit alpha
Tag	N-His
Concentration	>0.5 mg/ml
Purity	>95% by SDS-PAGE
Purification	AC
Protein Sequence	30-250/553
MW	28 kDa
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Expression System	E.coli
Endotoxins	Not analyzed
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)