

Product Datasheet

Human PLA2G7 Protein (orb1516755)

Description	Recombinant Human PLA2G7 Protein
Endotoxins	Not analyzed
Target	PLA2G7
Preservatives	8 M Urea (pH=9.0) with 1 mM DTT
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	22-441/441
Purity	> 90% as determined by SDS-PAGE
MW	48 kDa
Uniprot ID	Q13093
Expiration Date	6 months from date of receipt.