

Product Datasheet

Insulin Detemir Protein (orb1516882)

Description	Recombinant Insulin Detemir Protein
Endotoxins	Not analyzed
Target	169148-63-4
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	Tag free
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	Unknown
Purity	> 95% as determined by SDS-PAGE
MW	5.9 kDa
CAS Number	169148-63-4
Uniprot ID	P01308
Expiration Date	6 months from date of receipt.