

Product Datasheet

Human VEGF Protein (orb1516906)

Description	Recombinant Human VEGF Protein
Endotoxins	Not analyzed
Target	VEGFA
Form/Appearance	Liquid
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	27-232/232
Purity	> 90% as determined by SDS-PAGE
MW	25 kDa
Uniprot ID	P15692
Expiration Date	6 months from date of receipt.