

## Product Datasheet

### Human Insulin Protein (orb1516933)

<b>Description</b>	Recombinant Human Insulin Protein
<b>Endotoxins</b>	< 1.0 EU/μg as determined by endotoxin.
<b>Target</b>	INS
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	Tag free
<b>Note</b>	For research use only
<b>Protein Sequence</b>	1-110/110
<b>Purity</b>	> 95% as determined by SDS-PAGE
<b>MW</b>	6 kDa
<b>CAS Number</b>	11061-68-0
<b>Uniprot ID</b>	<b>P01308</b>
<b>Expiration Date</b>	6 months from date of receipt.