

## Product Datasheet

### Tacrolimus-BSA (orb1517031)

<b>Description</b>	Tacrolimus-BSA
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Coupling protein
<b>Hazard Information</b>	Non-Toxic
<b>Expiration Date</b>	12 months from date of receipt.