

Product Datasheet

SARS CoV2 Spike protein-RBD antibody (orb1517274)

Description	Goat polyclonal antibody to SARS CoV2 Spike protein-RBD
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SARS CoV2 Spike protein-RBD
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	0.2mg in 0.5mg/ml
Note	For research use only
Application notes	WB: Approx. 30 kDa band observed in recombinant RBD protein. Recommended concentration: 1µg/ml. Primary incubation was 1 hour.
Protein Sequence	10 aa long, RBD of spike protein
MW	30 kDa
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:2000.
Expiration Date	12 months from date of receipt.