

Product Datasheet

HA-Tag Antibody (orb1517280)

Description	HA-Tag Antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	Synthetic Peptide
Target	HA-Tag
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.