

Product Datasheet

V5-Tag Antibody (orb1517282)

| | |
|----------------------------|--|
| Description | V5-Tag Antibody |
| Species/Host | Mouse |
| Reactivity | All |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IF, IP, WB |
| Immunogen | Synthetic Peptide |
| Target | V5-Tag |
| Preservatives | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Monoclonal |
| Expiration Date | 12 months from date of receipt. |