

## Product Datasheet

### VSV-G-Tag Antibody (orb1517283)

<b>Description</b>	VSV-G-Tag Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	All
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IP, WB
<b>Immunogen</b>	Synthetic Peptide
<b>Target</b>	VSV-G-Tag
<b>Preservatives</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.