

Product Datasheet

African swine fever virus Phosphoprotein p30 Antibody (orb1517338)

Description	African swine fever virus Phosphoprotein p30(Ba71V-93) Monoclonal Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Phosphoprotein p30 protein (1-204AA)
Target	African swine fever virus Phosphoprotein p30
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P34204
Expiration Date	12 months from date of receipt.